

ON-Q[®] with Select-A-Flow[®]

P A T I E N T G U I D E L I N E I N S E R T

IMPORTANT: Start by reading the Patient Guideline folder then continue with this insert for complete information.

SELECT-A-FLOW RATE _____



SELECT-A-FLOW



The Select-A-Flow is a flow controller that allows your doctor to adjust the amount of medicine you receive to best meet your needs. Your Select-A-Flow may be either 1-7 or 2-14 ml/hr.

WARNING: Do not change the flow rate on the Select-A-Flow dial unless instructed by your doctor. Changing the flow rate without your doctor's instruction may result in the wrong dose of medicine delivered which could cause serious injury.

✓ **Do not tape tubing to skin.**
The Select-A-Flow should be worn outside your clothing.



24-HOUR PRODUCT SUPPORT

800.444.2728

VISIT US ON THE WEB

www.AskYourSurgeon.com